



Widener University

YOUR OFFICIAL WIDENER UNIVERSITY SUSSEX FINE QUALITY REGALIA MADE EASY

TAKE PRIDE IN YOUR ACADEMIC ACHIEVEMENT

By owning this specially designed doctoral regalia for Widener University, tailored to your exact specifications, your ensemble represents excellence in style and quality.



GOWN DETAILS:

- Revive sustainable light blue fabric gown
- Classic yoke with fluting on shoulders and back
- Decorative cord and button on back yoke
- Velvet chevrons and front panels
- Concealed center zipper
- Satin-lined bell sleeves
- Left pocket and access on right side
- Metallic gold piping on velvet chevrons and panels

HOOD DETAILS:

- Revive sustainable light blue fabric
- Black velvet
- Satin lining in school colors
- Cord and button
- Hood attachment cord

TAM DETAILS:

- Black velvet
- Deluxe 6-sided
- Gold bullion tassel

INCLUDED WITH PURCHASE:

- Initials on gown
- Garment bag and wooden hanger




REGALIA MEASURING

HOW TO MEASURE

Our gowns are designed in accordance with the Academic Costume Code and are custom tailored to fit.

ALL MEASUREMENTS ARE REQUIRED TO MAKE YOUR CUSTOM TAILORED REGALIA.



HEIGHT:

Record total height (in inches) from top of head to floor without shoes.

WEIGHT:

HEAD:

Enter your head circumference measurement (in inches). Measure 1" above the ears. Find and enter the corresponding cap or tam size from the chart below.

CAP/TAM SIZE	CAP SIZE	HEAD MEASUREMENT
XS	6 ¹ / ₄ " – 6 ¹ / ₂ "	19 ¹ / ₄ " – 20 ¹ / ₈ "
S	6 ⁵ / ₈ " – 6 ³ / ₄ "	20 ¹ / ₄ " – 21 ¹ / ₈ "
M	6 ⁷ / ₈ " – 7 ¹ / ₄ "	21 ¹ / ₄ " – 22 ⁷ / ₈ "
L	7 ³ / ₈ " – 7 ⁵ / ₈ "	23" – 24 ¹ / ₈ "
XL	7 ³ / ₄ " – 8 ¹ / ₈ "	24 ¹ / ₄ " – 26"

HIGH CHEST:

Enter your high chest measurement (in inches). Standing straight with your arms hanging straight by your sides, have a friend measure around your chest at your under arm. Tape measure should be straight around the back of your body.

FULL CHEST:

Enter your full chest measurement (in inches). Standing straight with your arms hanging straight by your sides, have a friend measure around your chest at the fullest point. Tape measure should be straight around the back of your body.

SHOULDER:

Enter your shoulder width (in inches) measuring across the back from tip of shoulder to tip of shoulder.

ARM:

Measure with your arms at your sides. Measure from the center back neck to shoulder bone, down the outside of the arm to the halfway point between your wrist and the first knuckle of your thumb.

NAME

EMAIL

PHONE

Email completed form to the Widener Bookstore at 1685mgr@follett.com with the subject line "Custom Doctoral Regalia Order - YOUR NAME." Upon receipt, a member of the bookstore staff will contact you to confirm your order. If you have questions, please call the bookstore directly at 610-490-7012.

Jostens[®]

JOSTENS.COM

©2021 Jostens, Inc. 212899